

Available Accessories for Models: MCO-17AC, MCO-34AC
Roller Base
Order Number: MCO-18RB
Description: For use in single or stacked installations. Solid steel base includes positioning plates for incubator levelers. High-impact casters permit easy location. Adjustable front mounting pins extend to floor to prevent movement when installation is complete. Pins retract if roller base must be moved.
Applies to Models: MCO-17AC, MCO-18AIC, MCO-18AIC(UV), MCO-18M, MCO-34AC, MCO-36AIC, MCO-36AIC(UV), MCO-36M
Shelf and Brackets
Order Number: MCO-46ST
Description: Additional shelves may be ordered. Each shelf is made of SANYO inCusaFe copper-enriched stainless steel. Includes two shelf brackets that install without tools.
Applies to Models: MCO-17AC, MCO-18AIC, MCO-18AIC(UV), MCO-18M, MCO-34AC, MCO-36AIC, MCO-36AIC(UV), MCO-36M
Stacking Kit
Order Number: MCO-17PS
Description: Stacking kit, for installing one incubator on top of another.
Applies to Models: MCO-17AC, MCO-34AC
Two Stage Regulator
Order Number: MCO-100L
Description: Two-stage gas regulator monitors cylinder supply and meters gas to incubator input. CGA fitting 320 for CO2 cylinder.
Applies to Models: MCO-17AC, MCO-18AIC, MCO-18AIC(UV), MCO-20AIC, MCO-34AC, MCO-36AIC, MCO-36AIC(UV), MCO-40AIC, MCO-5M