

## Available Accessories for Model: MCO-18M, MCO-36M

### Internal CO2 Cylinder Tank Switching Device

**Order Number:** MCO-21GC

**Description:**

**Applies to Models:**

MCO-18M, MCO-36M

### Optional UV Upgrade Kit

**Order Number:** MCO-18UVS2

**Description:**

**Applies to Models:**

MCO-18M, MCO-36M, MCO-5AC, MCO-5M

### Pad Lock Hasp

**Order Number:** MCO-18PK

**Description:**

**Applies to Models:**

MCO-18AIC, MCO-18AIC(UV), MCO-18M, MCO-36AIC, MCO-36AIC(UV), MCO-36M

### Roller Base

**Order Number:** MCO-18RB

**Description:** For use in single or stacked installations. Solid steel base includes positioning plates for incubator levelers. High-impact casters permit easy location. Adjustable front mounting pins extend to floor to prevent movement when installation is complete. Pins retract if roller base must be moved.

**Applies to Models:**

MCO-17AC, MCO-18AIC, MCO-18AIC(UV), MCO-18M, MCO-34AC, MCO-36AIC, MCO-36AIC(UV), MCO-36M

### Shelf and Brackets

**Order Number:** MCO-46ST

**Description:** Additional shelves may be ordered. Each shelf is made of SANYO inCusaFe copper-enriched stainless steel. Includes two shelf brackets that install without tools.

**Applies to Models:**

MCO-17AC, MCO-18AIC, MCO-18AIC(UV), MCO-18M, MCO-34AC, MCO-36AIC, MCO-36AIC(UV), MCO-36M