

Available Accessories for Model: MCO-20AIC

Communications Port

Order Number: MTR-480

Description: Located at rear of chamber, RS232/RS485 data port acquires information from microprocessor controller including temperature, CO2, door ajar status and humidity pan water level. Connector, cable and software not supplied.

Applies to Models:

MCO-18AIC, MCO-18AIC(UV), MCO-20AIC, MCO-36AIC, MCO-36AIC(UV), MCO-40AIC

Cylinder Switchover System, CO2

Order Number: MCO-21GC

Description: Automatically changes from primary to secondary gas cylinder when first cylinder is depleted. Audible alarm and flashing indicator on main control panel notifies user when switch has occurred. Warehouse installed or field installed by authorized service personnel only.

Applies to Models:

MCO-18AIC, MCO-18AIC(UV), MCO-20AIC, MCO-36AIC, MCO-36AIC(UV), MCO-40AIC

Data Acquisition Software

Order Number: MTR-2000

Description: For monitoring and/or controlling microprocessor system. Windows based, for installation on PC. LAN compatible, configurable SMTP server to internet for user PC or mobile phone delivery.

Applies to Models:

MCO-18AIC, MCO-18AIC(UV), MCO-20AIC, MCO-36AIC, MCO-36AIC(UV), MCO-40AIC

Independent Inner Door Kit

Order Number: MCO-20ID

Description: High impact, clear plastic doors attach to interior inventory system behind glass inner door. Customer installed; directions included.

Applies to Models:

MCO-20AIC

Pad Lock Hasp

Order Number: MCO-20PK

Description:

Applies to Models:

MCO-20AIC

Roller Base

Order Number: MCO-20RB

Description: For use in single or stacked installations. Solid steel base includes positioning plates for incubator levelers. High-impact casters permit easy location. Adjustable front mounting pins extend to floor to prevent movement when installation is complete. Pins retract if roller base must be moved.

Applies to Models:

MCO-20AIC, MCO-40AIC

Shelf and Brackets

Order Number: MCO-58ST

Description: The MCO-20AIC incubator holds up to fifteen shelves. Five shelves are included with each incubator. Additional shelves may be ordered. Each shelf includes two shelf brackets which insert without tools.

Applies to Models:

MCO-20AIC, MCO-40AIC

Two Stage Regulator

Order Number: MCO-100L

Description: Two-stage gas regulator monitors cylinder supply and meters gas to incubator input. CGA fitting 320 for CO2 cylinder.

Applies to Models:

MCO-17AC, MCO-18AIC, MCO-18AIC(UV), MCO-20AIC, MCO-34AC, MCO-36AIC, MCO-36AIC(UV), MCO-40AIC, MCO-5M